

# Recombinant Human POU2F2

Catalog No: #GP10655



Package Size: #GP10655-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human POU2F2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-200 amino acids of human POU2F2
Target Name	POU class 2 homeobox 2
Other Names	OCT2, OTF2, Oct-2
Accession No.	Swissprot:P09086Gene Accession:BC006101
Uniprot	P09086
GeneID	5452;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a homeobox-containing transcription factor of the POU domain family. The encoded protein binds the octamer sequence 5'-ATTTGCAT-3', a common transcription factor binding site in immunoglobulin gene promoters. Several transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only