

Recombinant Human OMG

Catalog No: #GP10658



Package Size: #GP10658-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human OMG
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 124-418 amino acids of human oligodendrocyte myelin glycoprotein
Target Name	oligodendrocyte myelin glycoprotein
Other Names	OMGP
Accession No.	Swissprot:P23515Gene Accession:BC018050
Uniprot	P23515
GenID	4974;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is cell adhesion molecule contributing to the interactive process required for myelination in the central nervous system. Mostly expression in Oligodendrocytes and myelin tissues.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only