

Recombinant Human DNAJC3

Catalog No: #GP10666



Package Size: #GP10666-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DNAJC3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human DnaJ (Hsp40) homolog, subfamily C, member 3
Target Name	DnaJ (Hsp40) homolog, subfamily C, member 3
Other Names	P58, HP58, ERdj6, PRKRI, P58IPK
Accession No.	Swissprot:Q13217Gene Accession:BC047936
Uniprot	Q13217
GenID	5611;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein with multiple tetratricopeptide repeat (TPR) motifs as well as the highly conserved J domain found in DNAJ chaperone family members. It is a member of the tetratricopeptide repeat family of proteins and acts as an inhibitor of the interferon-induced, dsRNA-activated protein kinase (PKR).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only