

Recombinant Human PDCD4

Catalog No: #GP10687



Package Size: #GP10687-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PDCD4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 220-469 amino acids of human PDCD4
Target Name	programmed cell death 4
Other Names	H731
Accession No.	Swissprot:Q53EL6Gene Accession:BC026104
Uniprot	Q53EL6
GeneID	27250;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a tumor suppressor and encodes a protein that binds to the eukaryotic translation initiation factor 4A1 and inhibits its function by preventing RNA binding. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only