

Recombinant Human SERPINF1

Catalog No: #GP10693



Package Size: #GP10693-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SERPINF1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 243 amino acids of human serpin peptidase inhibitor, clade F (alpha-2 antiplasmin, pigment epithelium derived factor), member 1
Target Name	serpin peptidase inhibitor, clade F (alpha-2 antiplasmin, pigment epithelium derived factor), member 1
Other Names	OI6, OI12, PEDF, EPC-1, PIG35
Accession No.	Swissprot:P36955Gene Accession:BC000522
Uniprot	P36955
GenelD	5176;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the serpin family, although it does not display the serine protease inhibitory activity shown by many of the other serpin family members. The encoded protein is secreted and strongly inhibits angiogenesis. In addition, this protein is a neurotrophic factor involved in neuronal differentiation in retinoblastoma cells.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only