

## Recombinant Human PRPH

Catalog No: #GP10695



Package Size: #GP10695-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PRPH
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 100-405 amino acids of human PRPH
Target Name	peripherin
Other Names	NEF4; PRPH1
Accession No.	Swissprot:P41219Gene Accession:BC032703
Uniprot	P41219
GeneID	5630;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a cytoskeletal protein found in neurons of the peripheral nervous system. The encoded protein is a type III intermediate filament protein with homology to other cytoskeletal proteins such as desmin, and is a different protein than the peripherin found in photoreceptors. Mutations in this gene have been associated with susceptibility to amyotrophic lateral sclerosis.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only