

## Recombinant Human PHB

Catalog No: #GP10698



Package Size: #GP10698-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PHB
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	prohibitin
Other Names	PHB1
Accession No.	Swissprot:P35232Gene Accession:BC013401
Uniprot	P35232
GeneID	5245;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is evolutionarily conserved, and its product is proposed to play a role in human cellular senescence and tumor suppression. Antiproliferative activity is reported to be localized to the 3' UTR, which is proposed to function as a trans-acting regulatory RNA. Several pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only