

Recombinant Human RCHY1

Catalog No: #GP10700



Package Size: #GP10700-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RCHY1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-261 amino acids of human RCHY1
Target Name	ring finger and CHY zinc finger domain containing 1
Other Names	ZCHY; ARNIP; CHIMP; PIRH2; RNF199; ZNF363; PRO1996
Accession No.	Swissprot:Q96PM5Gene Accession:BC047393
Uniprot	Q96PM5
GeneID	25898;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene has ubiquitin ligase activity. It mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including tumor protein 53, histone deacetylase 1, and cyclin-dependent kinase inhibitor 1B, thus regulating their levels and cell cycle progression. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only