

# Recombinant Human PFN1

Catalog No: #GP10709



Package Size: #GP10709-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PFN1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-140 amino acids of human profilin 1
Target Name	Profilin 1
Other Names	ALS18
Accession No.	Swissprot:P07737Gene Accession:BC002475
Uniprot	P07737
GeneID	5216;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the profilin family of small actin-binding proteins. The encoded protein plays an important role in actin dynamics by regulating actin polymerization in response to extracellular signals. Deletion of this gene is associated with Miller-Dieker syndrome, and the encoded protein may also play a role in Huntington disease. Multiple pseudogenes of this gene are located on chromosome 1.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only