

Recombinant Human PFN4

Catalog No: #GP10710



Package Size: #GP10710-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PFN4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	profilin family, member 4
Accession No.	Swissprot:Q8NHR9Gene Accession:BC029523
Uniprot	Q8NHR9
GeneID	375189;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene belongs to the profilin family. It binds to actin and affects the structure of the cytoskeleton. Mostly expression in testes, pituitary gland, midbrain, bone marrow, vulva.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only