

Recombinant Human PRL

Catalog No: #GP10721



Package Size: #GP10721-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human PRL |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 29-227 amino acids of human prolactin |
| Target Name | prolactin |
| Accession No. | Swissprot:P01236Gene Accession:BC015850 |
| Uniprot | P01236 |
| GeneID | 5617; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only