

Recombinant Human RETN

Catalog No: #GP10744



Package Size: #GP10744-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RETN
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 19-108 amino acids of human resistin
Target Name	resistin
Other Names	ADSF, RSTN, XCP1, FIZZ3, RETN1
Accession No.	Swissprot:Q9HD89Gene Accession:BC069302
Uniprot	Q9HD89
GeneID	56729;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene belongs to the family defined by the mouse resistin-like genes. The characteristic feature of this family is the C-terminal stretch of 10 cys residues with identical spacing. The mouse homolog of this protein is secreted by adipocytes, and may be the hormone potentially linking obesity to type II diabetes. Alternatively spliced transcript variants encoding the same protein have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only