

## Recombinant Human SCG3

Catalog No: #GP10752



Package Size: #GP10752-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | Recombinant Human SCG3  |
| Brief Description     | Recombinant Protein   |
| Immunogen Description | Fusion protein corresponding to a region derived from 20-232 amino acids of human secretogranin III |
| Target Name           | Secretogranin III   |
| Other Names           | SGIII   |
| Accession No.         | Swissprot:Q8WXD2Gene Accession:BC014539   |
| Uniprot               | Q8WXD2  |
| GeneID                | 29106;  |
| Storage               | -20~-80°C, pH 7.6 PBS   |

## Background

The protein encoded by this gene is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. Granins may serve as precursors for biologically active peptides. Some granins have been shown to function as helper proteins in sorting and proteolytic processing of prohormones; however, the function of this protein is unknown. Two transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only