

# Recombinant Human SOCS5

Catalog No: #GP10775



Package Size: #GP10775-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SOCS5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human suppressor of cytokine signaling 5
Target Name	suppressor of cytokine signaling 5
Other Names	CIS6; CISH6; Cish5; SOCS-5
Accession No.	Swissprot:O75159Gene Accession:BC032862
Uniprot	O75159
GeneID	9655;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene contains a SH2 domain and a SOCS BOX domain. The protein thus belongs to the suppressor of cytokine signaling (SOCS) family, also known as STAT-induced STAT inhibitor (SSI) protein family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling. The specific function of this protein has not yet been determined. Two alternatively spliced transcript variants encoding an identical protein have been reported.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only