

Recombinant Human SOCS3

Catalog No: #GP10776



Package Size: #GP10776-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SOCS3
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	suppressor of cytokine signaling 3
Other Names	CIS3; SSI3; ATOD4; Cish3; SSI-3; SOCS-3
Accession No.	Swissprot:O14543Gene Accession:BC060858
Uniprot	O14543
GeneID	9021;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the STAT-induced STAT inhibitor (SSI), also known as suppressor of cytokine signaling (SOCS), family. SSI family members are cytokine-inducible negative regulators of cytokine signaling. The expression of this gene is induced by various cytokines, including IL6, IL10, and interferon (IFN)-gamma. The protein encoded by this gene can bind to JAK2 kinase, and inhibit the activity of JAK2 kinase.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only