

Recombinant Human RPAP3

Catalog No: #GP10782



Package Size: #GP10782-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RPAP3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 185 amino acids of human RNA polymerase II associated protein 3
Target Name	RNA polymerase II associated protein 3
Accession No.	Swissprot:Q9H6T3Gene Accession:BC056415
Uniprot	Q9H6T3
GeneID	79657;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes an RNA polymerase II-associated protein. The encoded protein may function in transcriptional regulation and may also regulate apoptosis. Alternatively spliced transcript variants have been described. Forms an interface between the RNA polymerase II enzyme and chaperone/scaffolding protein, suggesting that it is required to connect RNA polymerase II to regulators of protein complex formation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only