

Recombinant Human SPARC

Catalog No: #GP10785



Package Size: #GP10785-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human SPARC |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 18-303 amino acids of human secreted protein, acidic, cysteine-rich (osteonectin) |
| Target Name | secreted protein, acidic, cysteine-rich (osteonectin) |
| Other Names | ON |
| Accession No. | Swissprot:P09486Gene Accession:BC004974 |
| Uniprot | P09486 |
| GenID | 6678; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes a cysteine-rich acidic matrix-associated protein. The encoded protein is required for the collagen in bone to become calcified but is also involved in extracellular matrix synthesis and promotion of changes to cell shape. The gene product has been associated with tumor suppression but has also been correlated with metastasis based on changes to cell shape which can promote tumor cell invasion.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only