

# Recombinant Human SPIN4

Catalog No: #GP10787



Package Size: #GP10787-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SPIN4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	spindlin family, member 4
Accession No.	Swissprot:Q56A73Gene Accession:BC112360
Uniprot	Q56A73
GeneID	139886;
Storage	-20~-80°C, pH 7.6 PBS

## Background

SPIN4 (spindlin family, member 4), is a 249 amino acids exact function is not known.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only