

Recombinant Human SYT6

Catalog No: #GP10794

Package Size: #GP10794-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SYT6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 200 amino acids of human synaptotagmin VI
Target Name	synaptotagmin VI
Other Names	sytVI
Accession No.	Swissprot:Q5T7P8Gene Accession:BC044948
Uniprot	Q5T7P8
GeneID	148281;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the synaptotagmin family. Synaptotagmins share a common domain structure that includes a transmembrane domain and a cytoplasmic region composed of 2 C2 domains, and are involved in calcium-dependent exocytosis of synaptic vesicles. This protein has been shown to be a key component of the secretory machinery involved in acrosomal exocytosis. Alternatively spliced transcript variants have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only