

ITLN2 Antibody

Catalog No: #42927

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	ITLN2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ITLN2 protein.
Immunogen Description	Fusion protein of human ITLN2
Target Name	ITLN2
Other Names	HL2; HL-2
Accession No.	Swiss-Prot#: Q8WWU7 Gene ID: 142683
Uniprot	Q8WWU7
GeneID	142683;
Calculated MW	36kd
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Immunohistochemistry: 1:100-1:200

Images



Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: Mouse small intestines tissue


Primary antibody: 1/300 dilution

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



Immunohistochemical analysis of paraffin-embedded Human Colorectal cancer tissue using #42927 at dilution 1/200,



Immunohistochemical analysis of paraffin-embedded Human Esophagus cancer tissue using #42927 at dilution 1/200,

Background

May play a role in the defense system against pathogens. Expressed only in the small intestine.

Note: This product is for in vitro research use only