

Recombinant Human SPATA7

Catalog No: #GP10800



Package Size: #GP10800-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SPATA7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 350-599 amino acids of human SPATA7
Target Name	spermatogenesis associated 7
Other Names	HSD3; LCA3; HSD-3.1; HEL-S-296
Accession No.	Swissprot:Q9P0W8Gene Accession:BC008656
Uniprot	Q9P0W8
GeneID	55812;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene, originally isolated from testis, is also expressed in retina. Mutations in this gene are associated with Leber congenital amaurosis and juvenile retinitis pigmentosa. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only