

## Recombinant Human BIRC2

Catalog No: #GP10804



Package Size: #GP10804-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human BIRC2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human baculoviral IAP repeat containing 2
Target Name	baculoviral IAP repeat containing 2
Other Names	API1, MIHB, HIAP2, RNF48, cIAP1, Hiap-2, c-IAP1
Accession No.	Swissprot:Q13490Gene Accession:BC016174
Uniprot	Q13490
GeneID	329;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a member of a family of proteins that inhibits apoptosis by binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2, probably by interfering with activation of ICE-like proteases. This encoded protein inhibits apoptosis induced by serum deprivation and menadione, a potent inducer of free radicals. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only