

Recombinant Human SYK

Catalog No: #GP10805



Package Size: #GP10805-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human SYK |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to C terminal 250 amino acids of human spleen tyrosine kinase |
| Target Name | Spleen tyrosine kinase |
| Other Names | p72-Syk |
| Accession No. | Swissprot:P43405Gene Accession:BC011399 |
| Uniprot | P43405 |
| GeneID | 6850; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes a member of the family of non-receptor type Tyr protein kinases. This protein is widely expressed in hematopoietic cells and is involved in coupling activated immunoreceptors to downstream signaling events that mediate diverse cellular responses, including proliferation, differentiation, and phagocytosis. It is thought to be a modulator of epithelial cell growth and a potential tumour suppressor in human breast carcinomas. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only