

Recombinant Human SYVN1

Catalog No: #GP10806



Package Size: #GP10806-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SYVN1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human synovial apoptosis inhibitor 1, synoviolin
Target Name	Synovial apoptosis inhibitor 1, synoviolin
Other Names	DER3; HRD1
Accession No.	Swissprot:Q86TM6Gene Accession:BC030530
Uniprot	Q86TM6
GenID	84447;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein involved in endoplasmic reticulum (ER)-associated degradation. The encoded protein removes unfolded proteins, accumulated during ER stress, by retrograde transport to the cytosol from the ER. This protein also uses the ubiquitin-proteasome system for additional degradation of unfolded proteins. Sequence analysis identified two transcript variants that encode different isoforms.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only