

THEMIS2 Antibody

Catalog No: #42928

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	THEMIS2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total THEMIS2 protein.
Immunogen Description	Fusion protein of human THEMIS2
Target Name	THEMIS2
Other Names	ICB-1; C1orf38
Accession No.	Swiss-Prot#: Q5TEJ8Gene ID: 9473
Uniprot	Q5TEJ8
GeneID	9473;
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Esophagus tissue using #42928 at dilution 1/200,



Immunohistochemical analysis of paraffin-embedded Human Liver tissue using #42928 at dilution 1/200,

Background

May constitute a control point in macrophage inflammatory response, promoting LPS-induced TNF production.

Note: This product is for in vitro research use only