

OIT3 Antibody

Catalog No: #42931

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OIT3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total OIT3 protein.
Immunogen Description	Fusion protein of human OIT3
Target Name	OIT3
Other Names	LZP
Accession No.	Swiss-Prot#: Q8WWZ8 Gene ID: 170392
Uniprot	Q8WWZ8
GeneID	170392;
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Liver cancer tissue using #42931 at dilution 1/100.

Background

May be involved in hepatocellular function and development.

Note: This product is for in vitro research use only