

# Recombinant Human TRAF3IP3

Catalog No: #GP10847



Package Size: #GP10847-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TRAF3IP3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 240-436 amino acids of human TRAF3IP3
Target Name	TRAF3 interacting protein 3
Other Names	T3JAM
Accession No.	Swissprot:Q9Y228Gene Accession:BC130417
Uniprot	Q9Y228
GeneID	80342;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The gene encodes a protein that mediates cell growth by modulating the c-Jun N-terminal kinase signal transduction pathway. The encoded protein may also interact with a large multi-protein assembly containing the phosphatase 2A catalytic subunit. Alternative splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only