

# Recombinant Human STXBP2

Catalog No: #GP10880



Package Size: #GP10880-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human STXBP2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human syntaxin binding protein 2
Target Name	Syntaxin binding protein 2
Other Names	FHL5; UNC18B; Hunc18b; UNC18-2; pp10122; MUNC18-2
Accession No.	Swissprot:Q15833Gene Accession:BC002869
Uniprot	Q15833
GeneID	6813;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the STXBP/unc-18/SEC1 family. The encoded protein is involved in intracellular trafficking, control of SNARE (soluble NSF attachment protein receptor) complex assembly, and the release of cytotoxic granules by natural killer cells. Mutations in this gene are associated with familial hemophagocytic lymphohistiocytosis. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only