

Recombinant Human IL18RAP

Catalog No: #GP10882



Package Size: #GP10882-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IL18RAP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 378-599 amino acids of human interleukin 18 receptor accessory protein
Target Name	Interleukin 18 receptor accessory protein
Other Names	ACPL; CD218b; IL18RB; CDw218b
Accession No.	Swissprot:O95256Gene Accession:BC106765
Uniprot	O95256
GeneID	8807;
Storage	-20~-80°C, pH 7.6 PBS

Background

Interleukin 18 receptor accessory protein, also known as IL18RAP and CDw218b (cluster of differentiation 218b), is a human gene. The protein encoded by this gene is an accessory subunit of the heterodimeric receptor for IL18. This protein enhances the IL18 binding activity of IL18R1 (IL1RRP), a ligand binding subunit of IL18 receptor. The coexpression of IL18R1 and this protein is required for the activation of NF- κ B and MAPK8 (JNK) in response to IL18.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only