

Recombinant Human VWF

Catalog No: #GP10884



Package Size: #GP10884-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human VWF
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	von Willebrand factor
Other Names	VWD; F8VWF
Accession No.	Swissprot:P04275Gene Accession:BC022258
Uniprot	P04275
GeneID	7450;
Storage	-20~-80°C, pH 7.6 PBS

Background

The glycoprotein encoded by this gene functions as both an antihemophilic factor carrier and a platelet-vessel wall mediator in the blood coagulation system. It is crucial to the hemostasis process. Mutations in this gene or deficiencies in this protein result in von Willebrand's disease. An unprocessed pseudogene has been found on chromosome 22.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only