

# Recombinant Human SERPINA5

Catalog No: #GP10889



Package Size: #GP10889-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SERPINA5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 5
Target Name	Serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 5
Other Names	PCI; PAI3; PAI-3; PCI-B; PROCI; PLANH3
Accession No.	Swissprot:P05154Gene Accession:BC008915
Uniprot	P05154
GenID	5104;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a member of the serpin family of proteins, a group of proteins that inhibit serine proteases. This gene is one in a cluster of serpin genes located on the q arm of chromosome 14. This family member is a glycoprotein that can inhibit several serine proteases, including protein C and various plasminogen activators and kallikreins, and it thus plays diverse roles in hemostasis and thrombosis in multiple organs.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only