

## Recombinant Human ABCF3

Catalog No: #GP10909



Package Size: #GP10909-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ABCF3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 518-707 amino acids of human ATP-binding cassette, sub-family F (GCN20), member 3
Target Name	ATP-binding cassette, sub-family F (GCN20), member 3
Other Names	EST201864
Accession No.	Swissprot:Q9NUQ8Gene Accession:BC009253
Uniprot	Q9NUQ8
GenID	55324;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This protein may display an antiviral effect against flaviviruses such as west Nile virus (WNV) in the presence of OAS1B. ATP-binding cassette (ABC) transporters regulate the transport of a variety of physiologic substrates. But this protein lacks transmembrane domains and is probably not involved in transport.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only