

Recombinant Human ACOT7

Catalog No: #GP10930



Package Size: #GP10930-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ACOT7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-165 amino acids of human ACOT7
Target Name	acyl-CoA thioesterase 7
Other Names	ACT; ACH1; BACH; LACH; LACH1; hBACH; CTE-II
Accession No.	Swissprot:O00154Gene Accession:BC017365
Uniprot	O00154
GeneID	11332;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the acyl coenzyme family. The encoded protein hydrolyzes the CoA thioester of palmitoyl-CoA and other long-chain fatty acids. Decreased expression of this gene may be associated with mesial temporal lobe epilepsy. Alternatively spliced transcript variants encoding distinct isoforms with different subcellular locations have been characterized.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only