

Recombinant Human ACSBG1

Catalog No: #GP10935



Package Size: #GP10935-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ACSBG1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human acyl-CoA synthetase bubblegum family member 1
Target Name	acyl-CoA synthetase bubblegum family member 1
Other Names	BG; BG1; BGM; LPD; GR-LACS
Accession No.	Swissprot:Q96GR2Gene Accession:BC009289
Uniprot	Q96GR2
GeneID	23205;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene possesses long-chain acyl-CoA synthetase activity. It is thought to play a central role in brain very long-chain fatty acids metabolism and myelinogenesis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only