

Recombinant Human ACSM5

Catalog No: #GP10936



Package Size: #GP10936-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ACSM5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 208 amino acids of human ACSM5
Target Name	acyl-CoA synthetase medium-chain family member 5
Accession No.	Swissprot:Q6NUN0Gene Accession:BC016703
Uniprot	Q6NUN0
GeneID	54988;
Storage	-20~-80°C, pH 7.6 PBS

Background

The acyl-CoA synthetase medium-chain (ACSM) family is comprised of ACSM1, ACSM2A, ACSM2B, ACSM3, ACSM4 and ACSM5, which encode for enzymes catalyzing the activation of medium-chain length fatty acids. ACSM5 is a 579 amino acid protein has a broad substrate specificity and utilizes magnesium or manganese as a cofactor. ACSM5 is expressed in kidney and liver.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only