

Recombinant Human FHL5

Catalog No: #GP10937

Package Size: #GP10937-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FHL5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 32-284 amino acids of human four and a half LIM domains 5
Target Name	four and a half LIM domains 5
Other Names	ACT; dJ393D12.2; RP3-393D12.2
Accession No.	Swissprot:Q5TD97Gene Accession:BC021723
Uniprot	Q5TD97
GenID	9457;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is coordinately expressed with activator of cAMP-responsive element modulator (CREM). It is associated with CREM and confers a powerful transcriptional activation function. CREM acts as a transcription factor essential for the differentiation of spermatids into mature spermatozoa. There are multiple polyadenylation sites found in this gene. Alternative splicing results in multiple transcript variants encoding the same protein.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only