

Recombinant Human ACTL8

Catalog No: #GP10942



Package Size: #GP10942-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ACTL8
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human actin-like 8
Target Name	Actin-like 8
Other Names	CT57
Accession No.	Swissprot:Q9H568Gene Accession:BC028909
Uniprot	Q9H568
GeneID	81569;
Storage	-20~-80°C, pH 7.6 PBS

Background

This protein belongs to the actin family. Actin-like 8 is a protein in humans that is encoded by the ACTL8 gene. The high levels of expression in testis and pancreas and weak expression in placenta. The subcellular location of this protein is cytoplasm, cytoskeleton.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only