

Recombinant Human ACTR3B

Catalog No: #GP10945



Package Size: #GP10945-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ACTR3B
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 246 amino acids of human ARP3 actin-related protein 3 homolog B (yeast)
Target Name	ARP3 actin-related protein 3 homolog B (yeast)
Other Names	ARP11; ARP3BETA
Accession No.	Swissprot:Q9P1U1Gene Accession:BC015207
Uniprot	Q9P1U1
GenID	57180;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the actin-related proteins (ARP), which form multiprotein complexes and share 35-55% amino acid identity with conventional actin. The protein encoded by this gene may have a regulatory role in the actin cytoskeleton and induce cell-shape change and motility. Pseudogenes of this gene are located on chromosomes 2, 4, 10, 16, 22 and Y. Alternative splicing results in multiple transcript variants and protein isoforms.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only