

Recombinant Human ACTRT1

Catalog No: #GP10946



Package Size: #GP10946-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ACTRT1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human actin-related protein T1
Target Name	Actin-related protein T1
Other Names	AIP1; ARIP1; ARPT1; HSD27
Accession No.	Swissprot:Q8TDG2Gene Accession:BC014597
Uniprot	Q8TDG2
GeneID	139741;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein related to the cytoskeletal protein beta-actin. This protein is a major component of the calyx in the perinuclear theca of mammalian sperm heads, and it therefore likely functions in spermatid formation. This gene is intronless and is similar to a related gene located on chromosome 1. A related pseudogene has also been identified approximately 75 kb downstream of this gene on chromosome X.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only