

## Recombinant Human ACOT8

Catalog No: #GP10949



Package Size: #GP10949-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ACOT8
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	acyl-CoA thioesterase 8
Other Names	hTE; PTE1; PTE2; PTE-1; PTE-2; HNAACTE; hACTE-III
Accession No.	Swissprot:O14734Gene Accession:BC117155
Uniprot	O14734
GeneID	10005;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a peroxisomal thioesterase that appears to be involved more in the oxidation of fatty acids rather than in their formation. The encoded protein can bind to the human immunodeficiency virus-1 protein Nef, and mediate Nef-induced down-regulation of CD4 in T-cells.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only