

Recombinant Human ACYP1

Catalog No: #GP10950



Package Size: #GP10950-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ACYP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-99 amino acids of human ACYP1
Target Name	acylphosphatase 1
Other Names	ACYPE
Accession No.	Swissprot:P07311Gene Accession:BC035568
Uniprot	P07311
GeneID	97;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the acylphosphatase family. The encoded protein is a small cytosolic enzyme that catalyzes the hydrolysis of the carboxyl-phosphate bond of acylphosphates. Two isoenzymes have been isolated and described based on their tissue localization: erythrocyte (common) type acylphosphatase encoded by this gene, and muscle type acylphosphatase. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only