

Recombinant Human ADAM20

Catalog No: #GP10953



Package Size: #GP10953-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADAM20
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 211-400 amino acids of human ADAM metallopeptidase domain 20
Target Name	ADAM metallopeptidase domain 20
Accession No.	Swissprot:O43506Gene Accession:BC025378
Uniprot	O43506
GeneID	8748;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The expression of this gene is testis-specific.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only