

Recombinant Human ADAM2

Catalog No: #GP10956



Package Size: #GP10956-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADAM2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human ADAM metalloproteinase domain 2
Target Name	ADAM metalloproteinase domain 2
Other Names	CT15; FTNB; PH30; CRYN1; CRYN2; PH-30b
Accession No.	Swissprot:Q99965Gene Accession:BC064547
Uniprot	Q99965
GeneID	2515;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. This member is a subunit of an integral sperm membrane glycoprotein called fertilin, which plays an important role in sperm-egg interactions.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only