

Recombinant Human ADAM7

Catalog No: #GP10957



Package Size: #GP10957-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADAM7
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ADAM metalloproteinase domain 7
Other Names	EAPI; GP83; GP-83; ADAM 7; ADAM-7
Accession No.	Swissprot:Q9H2U9Gene Accession:BC058037
Uniprot	Q9H2U9
GeneID	8756;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the ADAMs family of zinc proteases. These transmembrane proteins play roles in multiple processes including cell signaling, adhesion and migration. The encoded protein lacks protease activity and may play roles in protein-protein interactions and cell adhesion processes including sperm-egg fusion. Mutations in this gene may be involved in the progression of melanoma.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only