

Recombinant Human ADAM17

Catalog No: #GP10958



Package Size: #GP10958-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADAM17
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ADAM metalloproteinase domain 17
Other Names	CSVP; TACE; NISBD; ADAM18; CD156B
Accession No.	Swissprot:P78536Gene Accession:BC062687
Uniprot	P78536
GeneID	6868;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biologic processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The protein encoded by this gene functions as a tumor necrosis factor-alpha converting enzyme; binds mitotic arrest deficient 2 protein; and also plays a prominent role in the activation of the Notch signaling pathway.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only