

Recombinant Human ADAM21

Catalog No: #GP10959



Package Size: #GP10959-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADAM21
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 409-681 amino acids of human ADAM21
Target Name	ADAM metallopeptidase domain 21
Other Names	ADAM31; ADAM 21
Accession No.	Swissprot:Q9UKJ8Gene Accession:BC109024
Uniprot	Q9UKJ8
GeneID	8747;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The expression of this gene expression is testis-specific.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only