

Recombinant Human ADAMTSL1

Catalog No: #GP10960



Package Size: #GP10960-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADAMTSL1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 116-329 amino acids of human ADAMTS-like 1
Target Name	ADAMTS-like 1
Other Names	C9orf94; PUNCTIN; ADAMTSR1; ADAMTSL-1
Accession No.	Swissprot:Q8N6G6Gene Accession:BC100788
Uniprot	Q8N6G6
GeneID	92949;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a secreted protein and member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motif) family. This protein lacks the metalloproteinase and disintegrin-like domains, which are typical of the ADAMTS family, but contains other ADAMTS domains, including the thrombospondin type 1 motif. This protein may have important functions in the extracellular matrix. Alternative splicing results in multiple transcript variants encoding distinct proteins.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only