

Recombinant Human ADAMDEC1

Catalog No: #GP10963



Package Size: #GP10963-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADAMDEC1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 206-340 amino acids of human ADAM-like, decysin 1
Target Name	ADAM-like, decysin 1
Other Names	M12.219
Accession No.	Swissprot:O15204Gene Accession:BC069582
Uniprot	O15204
GeneID	27299;
Storage	-20~-80°C, pH 7.6 PBS

Background

This encoded protein is thought to be a secreted protein belonging to the disintegrin metalloproteinase family. Its expression is upregulated during dendritic cells maturation. This protein may play an important role in dendritic cell function and their interactions with germinal center T cells.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only