

Recombinant Human ADCK4

Catalog No: #GP10965



Package Size: #GP10965-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADCK4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 192-423 amino acids of human aarF domain containing kinase 4
Target Name	aarF domain containing kinase 4
Accession No.	Swissprot:Q96D53Gene Accession:BC013114
Uniprot	Q96D53
GeneID	79934;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein with two copies of a domain found in protein kinases. The encoded protein has a complete protein kinase catalytic domain, and a truncated domain that contains only the active and binding sites of the protein kinase domain, however, it is not known whether the protein has any kinase activity. Multiple transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only