

# Recombinant Human ADIPOQ

Catalog No: #GP10970



Package Size: #GP10970-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ADIPOQ
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 19-244 amino acids of human ADIPOQ
Target Name	adiponectin, C1Q and collagen domain containing
Other Names	ACDC; ADPN; APM1; APM-1; GBP28; ACRP30; ADIPQTL1
Accession No.	Swissprot:Q15848Gene Accession:BC096308
Uniprot	Q15848
GeneID	9370;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene is expressed in adipose tissue exclusively. It encodes a protein with similarity to collagens X and VIII and complement factor C1q. The encoded protein circulates in the plasma and is involved with metabolic and hormonal processes. Mutations in this gene are associated with adiponectin deficiency. Multiple alternatively spliced variants, encoding the same protein, have been identified.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only